

News from Mercury Products

Tuesday, December 18, 2007



APD1000
3.2



APD1000
2.0



APD
1150



AirChair

Dear Reader,

There has been a rather large gap since our last Newsletter but we have been so busy with oxygen work and the development of the new AirChair, that time just seems to have evaporated.

However, we are now back on track and I would like to take the opportunity to wish all our readers the Compliments of the Season.

The APD1000 Pulse Dose system is now in use on 4 continents and judging by the number of refills we are doing, is showing great benefit for our customers. Versions are also in use at sea so the word is spreading – if you want 3 times the duration against conventional constant flow systems, and a breathing regime that leaves your users with oxygen that does not dry the throat and allows the taste buds to return, then the APD 1000 is for you.

Note that we have introduced the 2 litre version for regional aircraft and the 1150 for business aircraft.

The new AirChair is attracting a lot of interest and is enjoying its first orders. We are producing 2 versions – the 15" (380 mm) wide and 17" (432 mm) wide. The all up weight is 6 kgs and importantly, the backs folds down to the rear, allowing the carer to stand close to the user, enhancing a safe transfer to the aircraft seat. We are careful to call this a Transit chair for use on the aircraft, not a wheelchair. In order to make a lightweight chair, which will fit narrow aisles, compromises have to be made and we believe that the AirChair has achieved those.

If you would like further details please view our website www.mercuryproducts.co.uk or e.mail us on merprod@yahoo.co.uk